

AI-Enabled VOC Sensor Arrays for Infection Risk Stratification and Clinical Decision Support

Deyu Liu, Taejun Ko, Derese Getnet, Buyu Yeh, Jennifer Dootz, Feiqi Li, V. Sujith Sajja, Robert Hopkins, Xiaonao Liu
Tao Treasures LLC (d/b/a NanoBioFAB), Frederick, MD, USA

Introduction

Limitations in Infection Detection for Decision support

- Microbial culture technique is essential for organism identification and susceptibility testing:
 - Challenge:** Delayed or unavailable in austere and Role 1 and 2 care.
- Molecular panels are faster but target-limited set of microbes but can also detect dead microbes/pathogens
 - Challenge:** Requires specialized skillsets, power requirements and the physical footprint may not be feasible in austere environments.

Potential Solution for the Limitations of Unmet Clinical Need

- Rapid detection solution using rugged diagnostic sensors that can be rapidly tailored to theater-relevant pathogens.
- Provide earlier clinical decision support than current available diagnostics.
- An adjunct aid for medicine support in the large-scale combat environment of high-risk individuals.
- Validated AI-driven solutions for medics/clinicians to monitor, aid with decision support and potential triage support in far forward (role 1 and 2) environments.

Our Technology and Solution




Technology

- Patented Volatile Organic Compound sensing nanoarray technology, miniaturized, customizable to different environments and microbes.
- Nanosensory array is integrated into the iNose™ system
 - High-sensitivity and high throughput system with minimal footprint





Anticipated end-state design of the technology for decision support in role 1-2 environment

Advanced Prototype - iNose™




Nanosensors integrated into the sensor holder cover to detect VOCs on COTS 12-well plate

A multiplexed sensor array with up to nine channels, each incorporating a distinct functional sensing material




On the left, iNose™ is recording the VOCs from the infected samples that are obtained from clinics. On the right components of iNose™ platform

Experimental Workflow




Study Component	Materials / Conditions
Foundational pathogen panel	ESKAPEE species: <i>E. faecium</i> , <i>S. aureus</i> , <i>K. pneumoniae</i> , <i>A. baumannii</i> , <i>P. aeruginosa</i> , <i>E. coli</i> , and <i>E. cloacae</i> ; one <i>S. aureus</i> / <i>P. aeruginosa</i> mixed culture; sterile controls
Clinical-matrix model	Pooled, de-identified, uninfected human synovial fluid derived from 205 remnant specimens; enriched 1:1 and spiked with <i>S. aureus</i> , <i>S. epidermidis</i> , or <i>P. aeruginosa</i> at 1–10 ⁶ CFU/mL
Sensor hardware	Foundational 19-type sensing array; optimized six-sensor array; integrated incubation/headspace system; Compact 12-well fixture
Expanded internal program	~100 clinically relevant targets, including representative near-neighbor, mixed-culture, and selected fungal conditions

Results



Semi-Quantitative Burden Estimation via Onset defined Time-To-Detection (TTD)

- For all three organisms, TTD decreased as initial bacterial concentration increased.
- S. aureus* was detected at 2.3 hr at 10⁶ CFU/mL and 6.5 hr at 1 CFU/mL. Similar inverse relationships were observed for *S. epidermidis* and *P. aeruginosa*.
- This platform was extended to ~100 clinically relevant targets, including near-neighbor organisms, mixed bacterial cultures, and selected fungal pathogens.



Distinct VOC Fingerprints Supported Multiclass Pathogen Classification

- Seven ESKAPEE pathogens
- Data has 96% mean classification accuracy, and the mixed-culture class achieved 99.8% mean class-wise accuracy.

Conclusion | Impact to Warfighter

Conclusion

- The platform differentiated 17 clinically relevant pathogens, near-neighbor organisms, mixed cultures, and selected fungal pathogens.
- Data has 96% mean classification accuracy, and the mixed-culture class achieved 99.8% mean class-wise accuracy using low-power VOC signatures with <5% variance for key analytes.

Developmental Status of the Technology

- TRL 5/6 platform with portable prototype devices and preliminary evaluation on clinical synovial fluid samples
- Ongoing development is focused on actionable readouts within 12 hrs, broader pathogen coverage, and matrix expansion to urine and blood.

Impact to Warfighter:

- The customizable sensor-manufacturing and analytics workflow supports theater-specific deployment, expansion toward >100 pathogen targets
- Earlier decision support for wound infection surveillance, casualty prioritization, and empiric treatment selection in austere settings.

Applicability to Medical Roles of Care:

- iNose can support Role 1-4 infection screening and triage by providing rapid infection risk assessment and classification support
- Samples include Gram-positive vs Gram-negative signatures and in body fluids when culture is unavailable, delayed, or resource-intensive.

This project was funded by
DHA SBIR HT942525CE008
HT942525CE045

Disclaimer: The views, opinions and/or findings contained in this research/presentation/ publication are those of the author(s)/company and do not necessarily reflect the views of the Department of Defense and should not be construed as an official DoD/Army position, policy or decision unless so designated by other documentation. No official endorsement should be made. Reference herein to any specific commercial products, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the U.S. Government."

Contact us:

Dr. Xianao Liu, CEO
email: xiaonaoliu@nanobiofab.com

Or scan to reach us



© NanoBioFAB Proprietary

NANOBIOFAB