

Maximum Efficiency of Extracorporeal CO₂ Removal Requires Nine -to-One Sweep Gas Flow to Blood Flow Rates for Metabolic Intervention Benefits

Marina Batchinsky¹, C. Brooke Sons^{1,2}, Ivan P. Rubalcava^{1,2}, Jacob Callahan¹, Daniel Wendorff¹, Brendan M. Beely, RRT¹, George Harea, PhD¹, Teryn R. Roberts, PhD^{1,2}, Andriy I. Batchinsky, MD^{1,2}, Lawrence A. Renna, PhD^{1,2}

¹Autonomous Reanimation and Evacuation Research Institute and Innovation Center, The Geneva Foundation, San Antonio, TX, USA; ²Department of Translational Medicine, University of the Incarnate Word, School of Osteopathic Medicine, San Antonio, TX, USA

The opinions or assertions contained herein are the private views of the author and are not to be construed as official or as reflecting the views of U.S. Government, the Department of Defense, or The Geneva Foundation.

Introduction

- High sweep gas flow rates during low-flow extracorporeal carbon dioxide removal (ECCO₂R) therapy rapidly remove CO₂ and stabilize pH levels
- High “sweep gas” flow rates lead to high oxygenator membrane ventilation-perfusion (V/Q) ratios, yet optimal V/Q is not well established
- Membrane oxygenator “dead space” is the fraction of the membrane that is ventilated but not perfused, and gas exchange is impeded

Objective

- Investigate the relationship between high V/Q ratios and CO₂ removal efficiency in an ex vivo circulation model

Hypothesis

- CO₂ removal becomes inefficient at higher V/Q ratios

Methods

- Ex vivo circuit (Figure 1) consists of two sections: **1)** physiological simulator (Cardiohelp + HLS Membrane, Maquet, Rastatt, Germany) to model *in vivo* circulation, and metabolic activity; **2)** Gas-exchange membrane oxygenator (Novalung XLung, Fresenius Medical Care NA, Waltham, MA) with heat exchanger used for ECCO₂R
- Exhaust gas sampled and dehumidified (using Nafion sampling line); CO₂ measured using a nondispersive IR sensor (1Hz) (SprintIR®-W; Gas Sensing Solutions Ltd; Cumbernauld, UK)




Figure 1: Ex vivo circuit design with physiological simulator, ECLS membrane oxygenator, and tubing circulating 2 liters of donor swine blood. Blood was heated to 38 °C using external water heater and heat exchanger in oxygenator. Schematic created with Biorender.com

Methods

- Circuits filled with 2 L heparinized donor swine blood (Animal Technologies Inc., Tyler, TX; USDA EU 142/2011, ISO certified)
- During circulation: a blood flow and sweep gas flow rate target was set:
 - Blood flow targets: Q_b = 0.5, 1.5, and 2.0 L/min
 - O₂ sweep targets: V_g = 3, 10 and 14 L/min
- Target combinations set, allowed 10 min equilibration
- Blood pCO₂ (GEM4000; Instrumentation Laboratories) and exhaust gas CO₂ (PeCO₂) analysis performed; new targets selected, process repeats
- V_g / Q_b combinations with > 3 paired measurements were analyzed (Figure 2) resulting in V/Q ratios of 1.50, 6.00, 6.67, 9.33, 20.00, and 28.00; n = 54 paired measurements
- CO₂ exhaust and V/Q ratios were compared using linear regression and Pearson correlation coefficients (PCC)
- Optimal V/Q was assessed with Welch’s t-test; p < 0.05
- Oxygenator performance calculations:
 - Ventilation/dilution index $VDI = 1 - PeCO_2 / pre\ pCO_2$
 - Membrane dead space index $MDSI = 1 - PeCO_2 / post\ pCO_2$
 - Gas-side CO₂ clearance $\dot{V}CO_2 = PeCO_2 / [P_{atm} - PH_2O] \times V_g$

Results




Figure 2: Histograms detailing sample counts for exhaust CO₂, pre-membrane CO₂, sweep gas flow, and blood flow.

Results




Figure 3: PeCO₂ vs pre-membrane blood pCO₂ at various V/Q ratios.

- PeCO₂ scaled linearly with pre pCO₂
- Slope decreased from 0.50 to 0.06 as V/Q increased from 1.50 to 28.0
- PCC decreased as V/Q increased (see table)

V/Q	n	slope	p (PCC)
1.50	8	0.501	0.990*
6.00	11	0.304	0.963*
6.67	7	0.251	0.993*
9.33	8	0.198	0.652
20.00	5	0.057	0.805
28.00	8	0.066	0.086

* p<0.05




Figure 4: VDI vs V/Q

- VDI, as a function of V/Q, showed saturating behavior, with an asymptote approaching VDI~0.9
- Consistent with a mixed ventilation–dilution effect on the exhaust signal




Figure 5: MDSI vs V/Q

- The MDSI showed no significant correlation with V/Q




Figure 6: V̇CO₂ vs V/Q

- Peak V̇CO₂ occurred at V/Q=9.33 (V_g=14, Q_b=1.5, t=3.73, p=0.007)
- Decreasing or plateauing values at higher V/Q ratios

Conclusions

- Optimal V/Q ratio for CO₂ removal is 9.33; higher ratios provide limited benefit in this model
- CO₂ removal inefficiency is attributable to overventilation rather than inadequate perfusion

Acknowledgments

This research was supported, in part, by UIWSOM and sponsored by AREVA Research Institute.