



Scale-Up and Technology Transfer of Manufactured Red Blood Cell Production to a CDMO for Clinical Manufacturing Readiness

Dr Rob Thomas PhD, CTO, Safi Biotherapeutics
MHSRS, Biomanufacturing Resilience, Orlando, FL, Aug 3rd, 2026

Disclosures and Learning Objectives

Disclosures

- Dr. Rob Thomas has no relevant financial or non-financial relationships with ineligible companies to disclose relating to the content of this activity.
- The views expressed in this presentation are those of the author and do not necessarily reflect the official policy or position of the Department of War, nor the U.S. Government.
- Commercial support was not received for this activity.

Learning Objectives


- Explain how a stirred-tank manufacturing platform produces cultured red blood cells and maintains product fidelity across 1 L to 10 L scale-up and transfer to a US-based CDMO
- Compare the in vitro and in vivo functional characteristics of manufactured red cells with donor-derived red cells
- Differentiate the novel properties of manufactured red cells — including glycerol-free cryopreservation resilience — from conventional donor red cell products

Overview

- **The Manufacturing Platform for a Cultured Red Cell Product**
- **In Vitro Data Indicating Functionality and Quality of Manufactured Red Cells**
- **In Vivo (mouse) Data Indicating Functionality and Quality of Manufactured Red Cells**
- **Novel Properties of Manufactured Red Cells**
- **Future Opportunities**





Scaled manufacturing of mRBCs

- Bioreactor platform, currently operating at up to 10 L
- Scale developed for pre-clinical and dose-escalation supply (~1 unit/batch)
- Format supports control and further scalability
- Unit operations bioreactor to bags being finalized




Initial clinical scale-up achieved: mRBC production

Expanded from 1L to 10L bioreactor production systems with high process and product fidelity




Product	235a	CD71	CD45
R59_1L	97.1 %	1.1 %	0.1 %
R60_10L	98.1 %	4.6 %	0.2 %



*Process Development
clinical process enabling*








1L mRBCs



*Scaled for clinical
production*

**10L mRBCs 2025
1 mRBC unit**

Manufacture Ready: Scaled batch-to-batch consistency at US-based GMP-ready CDMO facility (1L, 10L)




- Site transfer completed with success in growth, yield and quality via tech transfer runs at the 1L and 10L bioreactor scale
- US-based supply chain and material GMP compatibility

In vitro analysis across scales and tech transfer


Red cell assays demonstrate specific identity and functionality key to clinical application

Identity Phenotype: Manufactured(m) vs Donor (d)




Phenotype Drug Product Release parameters:

- >97% CD235a positive; CD45 negative
- >90% Draq5 negative (enucleation)





Representative In Vitro Functional Analysis at Clinical Scale Manufacture

mRBC O₂ dissociation curves during storage: Normal Hb oxygen affinity demonstrated via p50 saturation point and calculated Hill coefficient: Target range: >2.4



mRBC Deformability: Provides evidence of capacity to pass microcapillaries (Target >0.4 EI at 30 Pa)



Functional CQA data generated at Functional Fluidics, Dr. Patrick Hines

Development of enhanced in vitro analysis

Example of enhanced specification tailored for manufactured cell attributes

Batch	Day	Elmax (def)	*SS½	**Ohyper	***Δ Def-Osmo	****n_50	*****P50 (μ)	Verdict
10L Reference Run (Safi)	0	0.477	1.75	463	0.026	3.31	31.7	Pass
10L Reference Run (Safi)	7	0.465	2.01	514	0.003	4.00	33.5	High p50
1L Reference Run (Safi)	0	0.467	1.63	446	0.022	2.30	22.4	Low n_50
1L Reference Run (Safi)	7	0.458	1.85	492	0.009	2.42	20.1	Pass
▲ 1L Negative Control Run	0	0.211	1.06	290	-0.169	1.86	26.4	Fail (multi-assay)
Donor 1	0	0.554	2.19	404	0.015	3.15	16.5	Reference
Donor 1	7	0.575	2.79	375	0.047	3.44	15.8	Reference
Donor 1	0	0.584	2.40	400	0.044	3.12	17.1	Reference
Donor 1	7	0.604	3.82	374	0.082	3.29	15.6	Reference
CDMO 1L (ARMI)	0	0.451	1.62	431	0.018	3.02		
CDMO 1L (ARMI)	7	0.450	1.59	496	-0.001	2.86		
CDMO 10L (ARMI)	0	In progress	In progress	In progress	In progress	In progress		
CDMO 10L (ARMI)	7	In progress	In progress	In progress	In progress	In progress		

*SS1/2 shows low shear deformability is better in manufactured cells (relative to Elmax) compared to donor cells


**Ohyper shows osmolality of max deformability; lower = more dehydration

***Delta of max deformability from repeated stress (Def) vs single pass (Osmo) identifies fragile cells

****Hill coefficient calculated to account for p50 heterogeneity in manufactured product; expect >2.4

*****p50 highly variable due to sensitivity to allosteric modifiers and their loss dependent on temperature


▲ Negative control generated by overexposure to stress inducing reagents during culture to show clear failure



Hill plot:
Exclude
MetHb then
distribution
of p50
(spread,
skew, or
bimodality) ×
intrinsic n ×
metabolic
ΔpO₂

*Omics-based understanding driving continual improvement


As mRBC mature (enucleation) they become donor like across key biological domains




- Indices created to functionally group metabolic attributes show expected trajectory of canonical red cell properties in maturing manufactured product towards donor reference
- Deconstruction of indices enables in depth comparison to donor product, mechanistic insight into storage related changes, and data driven process/product improvement

*Refrigerated mRBC in vivo circulation model

Equivalent circulatory half-life of donor vs 1L/10L *manufactured RBCs in Mouse Model




CD235a (native marker) FITC half-life (days)



n.s. No significant difference between donor sources ($p = 0.86$, CD235a; $p = 0.80$, Biotin PE) or between 1L and 10L manufacturing batches ($p > 0.10$ all pairwise)



Biotin (clinically compatible marker) PE half-life (days)



*Cells approx. 7-10 days stored at point of transfusion after transatlantic shipment.



mRBC robustness and homogeneity support glycerol-free cryopreservation

- mRBCs recover ~90% with preserved deformability versus <10% for donor cells under equivalent storage media conditions or mRBC in standard DMSO cryopreservatives
- mRBCs higher glutathione/ATP content relative to donor cells is a mechanistic basis for advantages



Cryopreserved mRBC in vivo circulation model

Equivalent circulatory half-life of cryopreserved manufactured RBCs in Mouse Model



Cryo = commercially available (DMSO free); Cryo Improved uses further Safi in house protective additives (DMSO free)

Group	CD235a HL (mean $\pm$ SD, days)	n (mice)
Fresh Manufactured	2.89 $\pm$ 0.72	15 (3 batches)
Donor	3.48 $\pm$ 0.32	10 (sourced locally at transfusion site)
Cryo	2.46 $\pm$ 0.77	5
Cryo Improved	2.86 $\pm$ 0.39	11 (2 batches)

Fresh Manufactured

Donor

Cryo

Cryo Improved

Opportunities and Future Development

Near Term Clinical

- This platform has the capability to deliver manufactured red cells at a clinically meaningful scale; partnerships (ARMIBioFAB, Sciperio Inc.) are ensuring regulatory readiness of the technology and facilitating transfer for GMP manufacture

Near Term Manufacture Development

- Metabolomics data indicates opportunities to tune upstream formulation for higher production; specific targets identified in metabolism, membrane composition and oxidative resilience
- Evaluation of re-engineered gas transfer to enhance Oxygen/Carbon Dioxide flux at higher upstream production intensity and increased scale
- Downstream and formulation are specific to manufactured red cells: modifications for both cryo and refrigerated show path to enhanced recovery and storability (eg membrane active agents, metabolic/redox support)

Summary

- **A stirred-tank platform now produces cultured red cells at approximately unit scale (10 L), with unit operations defined for clinical supply.**
- **Process fidelity is maintained across scale-up and transfer to a US-based, GMP-ready CDMO.**
- **Manufactured red cells demonstrate identity, function, and in vivo circulatory survival in mouse models supporting clinical safety and efficacy**
- **Manufactured cells show novel, tunable attributes — including glycerol-free cryopreservation resilience — supported by omics-based characterization.**
- **Metabolomics-guided development and manufacturing partnerships are advancing regulatory and manufacturing readiness toward the clinic.**

Acknowledgements

- Thank you for listening
- Collaborators: ARMI BioFabUSA Team (NH); Sciperio Inc. (FL); Angelo D'Alessandro (CU); J. Roback and M. Lancelot, Emory University; Dr Patrick Hines and Functional Fluidics Team
- Funders: MTEC, DoW

This work was supported in part by the Department of War, Defense Health Program Combat Casualty Care Research Portfolio through the Medical Technology Enterprise Consortium (MTEC) under Other Transaction Number HT9425-24-9-0040. All animal studies were conducted under ACURO protocol MT24001.016.e001. The views, opinions, and/or findings are those of the author(s) and should not be construed as an official Department of War, Defense Health Program position, policy, or decision unless so designated by other documentation.