

Noninvasive Detection of Tibial Stress Fractures Using Surface Acoustic Waves

MHSRS 2026 26-5788 *Integrating Holistic & Technological Approaches for Musculoskeletal Injury Prevention Session*

Corresponding Author: Anna Knight Ph.D. Anna.Knight@jhuapl.edu

Anna Knight, PhD¹, Connor Pyles MS¹, Eyal Bar Kochba, PhD¹, Nadeau Hahne, MS¹, Nikki Steiner, MS¹, Constantine Demetropoulos, PhD¹, Catherine Carneal MS¹

¹ Johns Hopkins Applied Physics Laboratory, Laurel, MD

Disclosure

- This work was funded by the Johns Hopkins University Applied Physics Laboratory (JHU/APL) Research and Exploratory Development Mission Area and Global Health Mission Area Independent Research and Development (IRAD) programs. The views expressed are those of the authors and do not necessarily reflect those of JHU/APL.
- A nonprovisional patent application covering aspects of this technology has been filed
- All authors have no conflicts of interest to report


Challenge: Tibial Stress Fractures (TSF) in Military

- **Stress fractures (SF), particularly in the tibia, are common to athletes and the military**
 - Affects 3 - 8% of recruits during training [1] due to significant increase in exercise activities with new shoes and gear
 - 86.5% of SFs occur in the tibia [1] resulting in a 21.1 week mean rehab [2]
 - Recruits are 4x more likely discharge after sustaining a stress fracture [3]
- **Currently, by the time TSF is detectable, it is already too late to mitigate impact on athletic performance**
 - Diagnosed by clinical x-ray with “dreaded black line” (~1 mm resolution)
 - Indication that irreversible injury already incurred

Earlier detection of TSF initiation, enabled by portable acoustic diagnostic technology, will help mitigate injuries and optimize their performance abilities





Stress fracture of tibia




<http://www.svuhradiology.ie/case-study/stress-fracture-of-tibia/>

Solution: Acoustic Bone Surface Waves for TSF Detection

- Portable non-invasive, non-ionizing surface acoustic waves can detect changes in bone structure with greater resolution than current x-ray systems





Simulations and surrogate testing used to select distances between transmit (Tx) and receive (Rx), and minimum detectable crack size: 2.25 MHz and 55 degree PDMS wedge allows sensing ~60 mm




Solution: Acoustic Bone Surface Waves for TSF Detection

- Portable non-invasive, non-ionizing surface acoustic waves can detect changes in bone structure with greater resolution than current x-ray systems





Simulations and surrogate testing used to select distances between transmit (Tx) and receive (Rx), and minimum detectable crack size: 2.25 MHz and 55 degree PDMS wedge allows sensing ~60 mm




Solution: Acoustic Bone Surface Waves for TSF Detection


- Portable non-invasive, non-ionizing surface acoustic waves can detect changes in bone structure with greater resolution than current x-ray systems




Features in received acoustic signals indicate bone structural changes




Simulations and surrogate testing used to select distances between transmit (Tx) and receive (Rx), and minimum detectable crack size: 2.25 MHz and 55 degree PDMS wedge allows sensing ~60 mm



Fracture Detection in Bone: Manual Cuts




- Increase in time of flight with increasing cut depth
- Sensitive to cuts as small as 0.3 mm




System sensitive to detection of fractures under 1 mm, competing with the resolution of planar radiography

Fracture Detection in Bone: 4 specimen to date





- Increase in time of flight with increasing cut depth
- Sensitive to cuts as small as 0.5 mm






Same slope/trend across 4 specimen to date

Fracture Detection in Bone: Fatigue Induced



- ~1,200 N load, 8 mm displacement
- 10,000 cycles at 1.5 Hz
- Acoustic data every 1,000 cycles
- Measured in intact human lower extremity


Fracture Detection in Bone: Fatigue Induced



- ~1,200 N load, 8 mm displacement
- 10,000 cycles at 1.5 Hz
- Acoustic data every 1,000 cycles
- Measured in intact human lower extremity


Coherence of surface acoustic wave degrades over time during biomechanically relevant fatigue loading

Fracture Detection in Bone: Fatigue Induced




- ~1,200 N load, 8 mm displacement
- 10,000 cycles at 1.5 Hz
- Acoustic data every 1,000 cycles
- Measured in intact human lower extremity

Coherence vs Number Cycles of Fatigue, by Channel



Coherence: Fatigue Over 10,000 Cycles




Coherence trend holds over two specimen to date




Fracture Detection in Bone: Fatigue Induced



MicroCT
Inspection of Bone
Structure





- ~1,200 N load, 8 mm displacement
- 10,000 cycles at 1.5 Hz
- Acoustic data every 1,000 cycles
- Measured in intact human lower extremity



Coherence metric is validated with microCT, meaning acoustic signature is a first measured non-invasive metric of bone fatigue

Developing a Deployable Solution



- SWAP: Previously hardware 30 lbs, Updated SWAP: ~3 lbs, runs on USB-C, processing works on a phone application
- Non-provisional patent on “Inv P07907: Tibial Fracture Detection with Surface Acoustic Waves”, filed June 9 2025 no. 19/232,568 with the U.S. Patent & Trademark Office

Next Steps


Research

- Optimize Sensing Technology
 - Improve SNR
 - Reduce SWAP
- Validate in Preclinical Models
 - Conduct additional PMHS studies (n = 4)
 - Establish baseline variability across tibia, subjects, operator
 - Correlate acoustic signal with crack size, validate with microCT


Operationalization

- Translate to Human Subjects
 - Establish clinical partnership
 - Conduct healthy volunteer studies
 - Longitudinal training cohort
- Develop Predictive Models of Acoustic signature mapping injury risk
- Enable Injury Prevention & Treatment
 - Inform individualized training
 - Guide intervention
 - Monitor healing & integrate therapeutic ultrasound

Actively seeking clinical, collaborative, and transition partners



Questions?




Anna Knight, Ph.D.
anna.knight@jhuapl.edu




JOHNS HOPKINS
APPLIED PHYSICS LABORATORY

Live Demo


Baseline



Uninjured



Injured




Live Demo


Baseline



Uninjured



Injured



Specimen Loading Protocol


Objective: Apply cyclic physiologic loads representative of running to human tibia specimen. Should be well-controlled to allow us to...

- Identify elastic wave signatures representing distal anteromedial cortical strain and damage accumulation
- Compare these signature across test specimen
- Avoid damaging equipment

Candidate Loading Protocol applied on per-specimen basis

1. Apply instrumentation (strain gauges, acoustic transducers)
2. Align and pot specimen in test fixture
3. Load once at 1 mm
4. Determine **reference displacement** of crosshead at which 1000 $\mu\epsilon$ anteromedial strain is reached
 - In one specimen this is 6 mm, in the other 8 mm
5. Sinusoidal cyclic loading using displacement control up to **reference displacement** at 1.5 Hz for up to 10,000 cycles
 - Pause every 1,000 cycles at max and min displacement to acquire acoustic data

Test Characteristics	Literature	APL Multiscale MSK-FEM
Comp. Force	2-12x BW ^{3,5,6,7,12}	8.6x BW
Bending Moment	10-100x BW-mm ⁵	51.2 BW-mm
Cortical Strain	500 – 5,000 $\mu\epsilon$ ^{1,2,3,4,11}	1,250 $\mu\epsilon$
Strain Rate	0.01 – 0.03 ϵ/s ^{4,8}	0.025 ϵ/s
Cycles to Stress Fracture	1k – 10k ^{2,3}	-





JOHNS HOPKINS
APPLIED PHYSICS LABORATORY